## 1511 SANIMED



DESCRIPTION	Concentrated surface disinfectant	
SPECIFICATIONS	Concentrated disinfectant cleaner has fungicidal, bactericidal and virucidal action, killing gram positive and gram negative bacteria, fungi, mould and viruses (encapsulated and non- encapsulated) in household and professional environments. Particularly suitable for disinfecting surfaces and equipment in the agro-food, hotel and catering sector. It is also suitable for busy areas like public bathrooms, clinics and hospital areas. Gets rid of 99.9% of bacteria, mould and viruses.	
HOW TO USE AND DOSING	Dilute SANIMED® in water and spray or wipe with a sponge or mop. Soak small objects to be disinfected. For disinfection of surfaces with possible direct food contact rinse with drinking water. For virucidal action apply product to cleaned surfaces in a 5% solution, or 250 ml per 5 litres of water. HOW TO DILUTE: For normal cleaning dilute SANIMED® in a 0.5% solution, or 25 ml per 5 litres of water. For disinfection dilute SANIMED® in a 1.5% solution, or 75 ml per 5 litres of water. Leave to work for at least 5 minutes for bactericidal action or at least 15 minutes for fungicidal and virucidal action.	H.A.C.C.P.
SPECIFICATION	Appearance: Transparent green liquid pH: 10.5 +/- 0.5 Specific weight [g/ml]: 1.000 +/- 0.020 100 g of product contains: Pure alkyldimethylbenzylammonium chloride 10 g, Coformulants (*), colourings and water q.s. to 100 g.	
	The data given above are typical manufacturing values and not a specification.	
ADDITIONAL INFORMATION	Microbiological activity test: EN1276 "Bactericidal activity with interfering substances", EN1650 "Fungicidal activity with interfering substances", EN13697 "Bactericidal/fungicidal activity surface test", EN14476 "Virucidal activity suspension test", EN16777 "Virucidal activity surface test".	
PRODUCT COMPATIBILITY	Following recommended use methods it can be used on all materials commonly found in kitchens/used to manufacture various utensils. If in doubt, test material compatibility on a small concealed part before use.	
USE & STORAGE PRECAUTIONS	The product is formulated for professional use only. Store in original containers. Avoid extreme temperatures. Consult safety sheet for further information on use and disposal.	
CERTIFICATIONS/ CONFORMITY	The surfactants contained in this product are biodegradable according to the criteria defined by Regulation (EC) 648/2004 and subsequent amendments.	
CHEMICAL COMPOSITION	Consult safety sheet for further information on product ingredients.	EDICO.











	NO	CODE	1511	EAN	8032680391248		
	LETIZATI	Pcs x Box	2	Volume L	5	Fl. Oz.	169,07
		Box x Layer	18	Peso Kg	5	Oz.	176,37
	PAL	Box x Pit	72				

## WWW.SANITECITALIA.COM

LAST UPDATE: 15/02/2021 REV. 0/2021
TALCHIMICA SRL - RIVIERA DEI MAESTRI DEL LAVORO, 10 - 35127 PADOVA (PD) - ITALY / PHONE +39 049 8792456 - FAX +39 049 8791424 - INFO@ITALCHIMICA.IT - WWWITALCHIMICA.IT