

# 2049HC

## QUANTUM OXYGEN



**SANITEC**  
IGIENE SICURA

<b>DESCRIPTION</b>	Concentrated active oxygen-based fabric bleach
<b>SPECIFICATIONS</b>	Powerful stain remover for colours also For use on all types of fabrics For use even at temperatures over 60°C Optimised for use with Quantum System products
<b>HOW TO USE AND DOSING</b>	For use exclusively with L-Brain Connect automatic dosing systems Must be used in the pre-wash or washing phase Recommended dose: 2-7 ml/kg of dry washing. The recommended dose refers to product use in optimal conditions. Consult with Sanitec Service technical support for advice on dosing adjustment.
<b>SPECIFICATION</b>	Appearance: transparent colourless liquid pH: 2.5 +/- 0.5 Specific weight [g/ml]: 1.130 +/- 0.020  The data given above are typical manufacturing values and not a specification.
<b>PRODUCT COMPATIBILITY</b>	Do not use on nylon garments. Optimal pH for bleach is 9.0-10.5. Good results are also obtained at 60°C and pH 11.0-11.5
<b>USE &amp; STORAGE PRECAUTIONS</b>	The product is formulated for professional use only. Store in original containers. Avoid extreme temperatures and exposure to sunlight. Consult safety sheet for further information on use and disposal.
<b>CERTIFICATIONS/ CONFORMITY</b>	The product does not contain surfactants
<b>CHEMICAL COMPOSITION</b>	Consult safety sheet for further information on product ingredients.



**PROFESSIONAL  
ACADEMY**



Discover further  
informations on your  
smartphone.



PALLETIZATION	<b>CODE</b>	2049HC	<b>EAN</b>	8054633835777		
	<b>Pcs x Box</b>	1	<b>Volume L</b>	15	<b>Fl. Oz.</b>	507,21
	<b>Box x Layer</b>	13	<b>Peso Kg</b>	16,95	<b>Oz.</b>	597,89
	<b>Box x Pit</b>	52				

**WWW.SANITECITALIA.COM**

LAST UPDATE: 15/02/2021

REV.00/2021