


2050HC QUANTUM OXYGEN



SANITEC
IGIENE SICURA

DESCRIPTION	Concentrated active oxygen-based fabric bleach
SPECIFICATIONS	Powerful stain remover for colours also For use on all types of fabrics For use even at temperatures over 60°C Optimised for use with Quantum System products
HOW TO USE AND DOSING	For use exclusively with L-Brain Connect automatic dosing systems Must be used in the pre-wash or washing phase Recommended dose: 2-7 ml/kg of dry washing. The recommended dose refers to product use in optimal conditions. Consult with Sanitec Service technical support for advice on dosing adjustment.
SPECIFICATION	Appearance: transparent colourless liquid pH: 2.5 +/- 0.5 Specific weight [g/ml]: 1.130 +/- 0.020 The data given above are typical manufacturing values and not a specification.
PRODUCT COMPATIBILITY	Do not use on nylon garments. Optimal pH for bleach is 9.0-10.5. Good results are also obtained at 60°C and pH 11.0-11.5
USE & STORAGE PRECAUTIONS	The product is formulated for professional use only. Store in original containers. Avoid extreme temperatures and exposure to sunlight. Consult safety sheet for further information on use and disposal.
CERTIFICATIONS/ CONFORMITY	The product does not contain surfactants
CHEMICAL COMPOSITION	Consult safety sheet for further information on product ingredients.



**PROFESSIONAL
QUALITY**



Discover further
informations on your
smartphone.



PALLETIZATION	CODE	2050HC	EAN	8032680390173		
	Pcs x Box	2	Volume L	5	Fl. Oz.	169,07
	Box x Layer	18	Peso Kg	5,65	Oz.	199,30
	Box x Pit	72				

WWW.SANITECITALIA.COM

LAST UPDATE: 15/02/2021

REV.00/2021