

2053HC QUANTUM CHLOR



SANITEC
IGIENE SICURA

DESCRIPTION	Concentrated active chlorine-based fabric bleach
SPECIFICATIONS	Powerful stain remover on most fabrics Effective at low temperatures Guarantees excellent cleaning Optimised for use with Quantum System products
HOW TO USE AND DOSING	For use exclusively with L-Brain Connect automatic dosing systems Must be used during pre-wash phase Recommended dose: 3-9 ml/kg of dry washing. The recommended dose refers to product use in optimal conditions. Consult Sanitec Service technical support for dosing adjustment
SPECIFICATION	Appearance: transparent pale yellow liquid pH: 13.0 +/- 0.5 Specific weight [g/ml]: 1.250 +/- 0.020 The data given above are typical manufacturing values and not a specification.
PRODUCT COMPATIBILITY	Do not use on nylon garments. Optimal pH for bleach is 9.5 at ~30°C. At pH 10.5 you can wash up to 50°C. Do not bleach chlorhexidine stains which may set permanently in the fibres. Completely neutralise bleach in the rinse phase to avoid yellowing and damage to garments during ironing.
USE & STORAGE PRECAUTIONS	The product is formulated for professional use only. Store in original containers. Avoid extreme temperatures and exposure to sunlight. Consult safety sheet for further information on use and disposal.
CERTIFICATIONS/ CONFORMITY	The product does not contain surfactants
CHEMICAL COMPOSITION	Consult safety sheet for further information on product ingredients.



**PROFESSIONAL
QUALITY**



Discover further
informations on your
smartphone.



PALLETIZATION	CODE	2053HC	EAN	8032680390180		
	Pcs x Box	2	Volume L	5	Fl. Oz.	169,07
	Box x Layer	18	Peso Kg	6,25	Oz.	220,46
	Box x Pit	72				

WWW.SANITECITALIA.COM

LAST UPDATE: 15/02/2021

REV.00/2021